

Intracellular Solution Protocol

3/14/2019

- Calibrate the microsometer
 - Turn on
 - Run diagnostic
 - Run calibration
 - Use calibration standard solution
 - 50 x3
 - 850 x3
 - 2000 x3
- Make stocks fresh **EVERYTIME**
 - Best way to measure stocks is to put the entire tube in to the weigher, then add ddH2O
- For pH
 - Calibrate the pH before use here every time
 - For **EPSC** and **IPSC** uses **cesium-hydroxide** to increase pH
 - For **K-ICS** use **potassium-hydroxide** to increase pH
 - Use **HCL** to lower pH for **all**
 - **pH must be 7.35**
- Keep all solutions below 50 mL, around 40 mL for pH and osmolality adjustment
- For osmole test you can only lower the osmoles so **CAREFULLY APPROACH AIM (290 mOSM)**
 - To adjust
 - $Current\ Osmole\ (363) \times current\ volume\ (43) = 15,609$
 - $\frac{15,609}{Aim\ Osmoles\ (290)} = 53.82$
 - $53.82 - current\ volume\ (43) = 10.824$
 - Add 10.824 μ L **BUT** add less than that as you can only go down with osmoles, once past that's it

EPSC	mM	MW	g/50ml	stock	Check	Check	Check
Cs-methanesulfonate (CsCH3SO3)	125	228	1.425				
CsCl	5	168.36	0.04209	2.5 ml			
HEPES	10	238.31	0.119155	0.5 ml			
phosphocreatine	10	453.4	0.2267				
EGTA	0.2	380.35	0.0038035	0.1 ml			
ATP Na	4	551.14	0.110228				
GTP Na	0.3	523	0.007845	0.75 ml			
QX314	5	343.3	0.085825				
MgCl2.6H2O	4	203.3	0.04066	0.2 ml			
IPSC	mM	MW	g/50ml				
Cs-methanesulfonate (CsCH3SO3)	90	228	1.026				
CsCl	40	168.36	0.33672				
HEPES	10	238.31	0.119155	0.5 ml			
phosphocreatine	10	453.4	0.2267				
EGTA	0.2	380.35	0.0038035	0.1 ml			
ATP Na	4	551.14	0.110228				
GTP Na	0.3	523	0.007845	0.75 ml			
QX314	5	343.3	0.085825				
MgCl2.6H2O	4	203.3	0.04066	0.2 ml			
K ICS	mM	MW	g/50ml				
K-gluconate	125	234.2	1.46375				
KCl	5	74.55	0.0186375	0.25 ml			
HEPES	10	238.31	0.119155	0.5 ml			
phosphocreatine Na2	10	453.4	0.2267				
EGTA	0.2	380.4	0.003804	0.1 ml			
ATP-Na 2	4	551.4	0.11028				
GTP-Na 3	0.3	523.2	0.007848	0.75 ml			
MgCl 2 •6H 2 O	4	203.3	0.04066	0.2 ml			
pH7.35, 290mOsm							
Stock (M)			g/10ml	Add in (ml)			
0.1 CsCl			0.17	2.50			
1 HEPES			2.38	0.50			
1 Mgcl2			2.03	0.20			
0.1 EGTA			0.38	0.10			
0.02 GTP-Na			0.10	0.75			
1 KCl			0.7455	0.25			